



Discreet design.

User-friendly handle for good grip and non-touch technique.

Adhesive tab on back for reliable hanging.

Instant activation - sterile salt solution included, squeeze and it's ready to use.

Peel-open tab.

The outer packaging doubles as a hygienic and discreet disposal pouch after use.

Not manufactured with PVC, phthalates or latex.

Smooth eyelets.^{1,2}



LoFric catheters are coated with the unique Urotonic™ Surface Technology. This ensures the catheter remains smooth and safe, both during insertion and withdrawal.^{3,4}

Product description

LoFric® Sense™, single-use, hydrophilic urinary catheter with salt solution for activation, intended for intermittent urinary catheterisation.

Sterilisation and shelf life

Sterilised by irradiation. Use before expiry date on package.

Package

30 catheters per customer box. 4 customer boxes per shipping box.

Material composition

Hydrophilic surface coating: Polyvinyl pyrrolidone (PVP).

Catheter: Polyolefin-based elastomer (POBE).

Sachet: Laminate foil of polyethylene terephthalate (PET), aluminum and polyethylene (PE). Sterile water with sodium chloride (NaCl).

Item packaging: Laminate foil of polyethylene terephthalate (PET) and polyethylene (PE).

Product range

Reference	Type	Length	Size
416 08 25	LoFric® Sense™ Nelaton (female)	15 cm	CH 08
416 10 25	LoFric® Sense™ Nelaton (female)	15 cm	CH 10
416 12 25	LoFric® Sense™ Nelaton (female)	15 cm	CH 12
416 14 25	LoFric® Sense™ Nelaton (female)	15 cm	CH 14

The effective length is 10 cm. Effective length, length of the catheter that can be inserted into the body.



User-friendly packaging.



LoFric® Sense™

References

1. Lundgren et al., Spinal Cord, 2000
2. Newman et al., Springer, 2018
3. Waller et al., Spinal Cord, 1997
4. Dunger et al., Acta Endocrinologica, 1988

At Wellspect we develop innovative continence care solutions that improve quality of life for people with bladder and bowel problems. We inspire our users to build self-confidence and independence as well as good health and well-being. We have been leading the industry for over 40 years with our product brands LoFric® and Navina™. We always aim to minimise the environmental impact of our products and passionately strive to become climate neutral. We work together with users and healthcare professionals to improve clinical outcomes in a sustainable way, now and for the future.

Wellspect. A Real Difference.

For more information about our products and services, please visit [Wellspect.co.uk](https://www.wellspect.co.uk).

Join the conversation on LinkedIn, Twitter, Facebook and Instagram.

[wellspect.co.uk](https://www.wellspect.co.uk)

Wellspect Limited, Stroudwater Business Park, Brunel Way, Stonehouse, Gloucestershire, GL10 3GB, Great Britain.
Phone: 0800 652 3350. www.wellspect.co.uk

LoFric



Manufacturer: Wellspect HealthCare, Aminogatan 1, P.O. Box 14, SE-431 21 Mölndal, Sweden.
Phone: +46 31 376 40 00. www.wellspect.com

CE 2797

Rx Only